

Observation and Result:-

Test Type	Ppt appear	Time Taken
Seliwanoff's Test	Cherry red pink colored in 30 sec	30 seconds

Result:-

- Hexoketoses are confirmed
- So, Seliwanoff's test is a ketoses identifying test.

EXPERIMENT # 05

30-Nov-2020

Title :-

Determine the Seliwanoff's Test :-
(For keto Sugars)

Principle :-

This test can differentiate between aldose and ketose because aldohexoses will produce light pink color after sometime.

However after prolong boiling, the color is deepen, so restrict the prolong boiling.


This test is based on the fact that when carbohydrates are heated, ketoses are more rapidly dehydrated than aldoses.

Therefore boiling should before 30 seconds. In concentrated HCl, ketoses undergo dehydration to fulfil or to yield furfural derivatives more rapidly than do aldoses.

These derivatives forms complexes with resorcinol to yield deep cherry red color. It is tinned color reaction specific for ketoses.

Ketoses react more readily with Seliwanoff's reagent because: ketones have 3 carbon atoms where aldehyde has 2 hydrogen atoms and one carbon atoms.

Formula :-



Fructose + HCl \rightarrow Hydroxymethyl Furfural + Levulinic Acid
 Hydroxymethyl Furfural + Resorcinol \rightarrow Cherry red complex
 Sucrose will also give's a positive test

Preparation of Seliwanoff's Reagent :-

Dissolve 0.05 grams of resorcinol in 50ml of distilled water. Add 33ml of Conc. HCl (37%) to this solution and further add distilled water upto total volume 100ml.

Procedure:-

Take 3ml of seliwanoff's reagents and one ml of given carbohydrate solution in a test tube and mix.

- Boil for 30 seconds only and then cool the solution. Note the color that appears.
- Appearance of a cherry red color or pink color within 30 sec indicates the presence of ketohexoses.

Precautions :-

1. The final concentration of HCl in Seliwanoff's should not exceed 12% HCl, because in the presence of very strong acid aldolases may be

converted to ketohexoses and give a false positive test.

2. Similarly in the case with the boiling and hence should not be done for a prolonged period. It goes without saying that other precautionary measures of the laboratory have to be taken.

Clinical Correlation:-

- HCL Poisoning
 - Carbon monoxide poisoning
-